(AU)
OUR GOODS COME WITH GUARANTEES THAT CANNOT BE EXCLUDED UNDER THE AUSTRALIAN CONSUMER LAW. YOU ARE ENTITLED TO A REPLACEMENT OR REFUND FOR A MAJOR FAILURE AND COMPENSATION FOR ANY OTHER REASONABLY FORSEEABLE LOSS OR DAMAGE. YOU ARE ALSO ENTITLED TO HAVE THE GOODS REPAIRED OR REPLACED IF THE GOODS FAIL TO BE OF ACCEPTABLE QUALITY AND THE FAILURE DOES NOT AMOUNT TO A MAJOR FAILURE.

US

One (1) Year Limited Consumer Product Warranty

STYLIDEAS LTD ("STYLIDEAS") warrants that the product shall be free from defects in material and workmanship for a period of 1 year from the date of original purchase ("Limited Warranty Period"), provided the product was used in accordance with its use and care instructions (e.g., in the intended environment and under normal circumstances).

Your Limited Warranty Benefits

During the Limited Warranty Period and subject to this limited Warranty Policy, STYLIDEAS will, in its sole and exclusive discretion, either (i) refund the purchase price if the purchase was made directly from the online STYLIDEAS store, (ii) repair any defects in material or workmanship, (iii) replace the product with another product of equal or greater value, or (iv) provide store credit in the amount of the purchase price.

Claiming Your Limited Warranty Service in 5 Simple Steps:

1. Make sure your product is within the specified limited warranty period.
2. Make sure you have a copy of the invoice and order ID or proof-of-purchase.
3. Make sure you have your product. DO NOT dispose of your product before contacting us.
4. Contact our Customer Support Team via +44(0)20 3225 5559 teamtom@styleideas.com
5. Once our Customer Support Team has approved your request, please return the product with a copy of the invoice and order ID.

©2025 STYLIDEAS Ltd. All rights reserved. STYLIDEAS owns trademarks, registered trademarks, patents, patents pending, design registrations and protection of the copyright of its products. www.styleideas.com/vp. Distributed by STYLIDEAS Ltd. UK: Amshold House, 1010 20W, UK. +44(0)20 3225 5559 teamtom@styleideas.com. EU: OBELIS S.A. Btl. General Wahnis 53, B-1030 Brussels, Belgium AUS: 6/7 Smeaton Grange Rd, Smeaton Grange, NSW, USA: 2020 Ocean Ave Ronkonkoma NY, 11779

NS01A

Made in China.

For indoor use only

STYLPRO™
Neck & Face Smoother

(EN) USE & CARE GUIDE

FOR YOUR SAFETY AND CONTINUED ENJOYMENT OF THIS PRODUCT, ALWAYS READ INSTRUCTION LEAFLET CAREFULLY BEFORE USING.

The product is intended for private use only and is not intended for commercial use or for use in other applications. The product is designed for the neck and face area, do not use the product on other parts of the body.

SPECIFICATION:


Name: STYLPRO™ Neck & Face Smoother
Model: NS01A
Device Size: 89mm(L) x 47mm(W) x 120mm(H)
Materials: Coated ABS Plastic
Nominal Voltage: 3.7V
Rated capacity of the battery: 750mAh
Energy Capacity of the battery: 2.77 Wh
Charge Voltage: 5V
Input power supply: USB socket, 5V, 1A
Charging Time: 2hrs

INCLUDED CONTENTS:

- 1 x STYLPRO™ Face and Neck Smoother
- 1 x Protective Cap
- 1 x Charging Cable
- 1 x Instruction Booklet

PARTS IDENTIFICATION:

1. Massage Head 2. Display 3. Power/Mode Button 4. Red/Green/Blue Light Indicator 5. battery charging indicator 6. USB-C charging port 7. Protective Cap 8. USB-C charging Cable.

SAFETY GUIDELINES:

1. We do not advise using this device if you are epileptic, photo-phobic, or light sensitive, if you are pregnant or breastfeeding or have had skin cancer in the last 12 months, have undergone facial plastic surgery unless under medical supervision, have had a skeletal reconstruction or titanium implant, use with care during your menstrual cycle.
2. For adult use only. Not intended for use by children. This is not a toy.
3. Wipe gently with a damp cloth.
4. Avoid storing the STYLPRO™ Fabulous Firmer in hot, humid, or damp areas.
5. If the STYLPRO™ Fabulous Firmer malfunctions or is damaged in any way, stop using it immediately. Contact STYLPRO™ CUSTOMER SERVICE (See section LIMITED TWO-YEAR WARRANTY for contact details).
6. Use LED lights with care. Extended use and prolonged exposure of LED lights WITHOUT EYE PROTECTION may cause damage. LEDs are not replaceable.
7. Do not use on sensitive, irritated or broken skin. Do not use on swollen or sunburnt skin or if you have any other skin conditions.
8. Do not look directly into the LED lights.

Test yourself for light sensitivity ON YOUR ARM prior to using THE DEVICE on your face:

1. Press the On/Off/SETTING button.
2. Place the device on your arm and conduct a 10 minute SESSION.
3. After 10 minutes, turn the device off. Visually check the area. THE AREA MAY INITIALLY TURN RED DUE TO THE HEATING OF YOUR SKIN BUT IF THE REDNESS LASTS FOR MORE THAN 2 hours, your skin has a sensitivity to light. SEEK MEDICAL ADVICE BEFORE DECIDING HOW BEST TO USE THIS DEVICE ON YOUR SKIN.

CHARGING:

Insert the USB-C cable (provided) into the port at the bottom of the Device. Attach to a USB-C friendly power source. The battery indicator light will be flash red when charging and stop flashing when fully charged.

CLEANING:

1. Please wipe off any excess skin care products which may accumulate in the massage head.
2. Wipe clean with a damp cloth only.
3. Do not submerge in water.
4. Please don't use cleaning product to clean.
5. Please store it in a cool dry place.
6. Please don't use if it is damaged.

DANGER: TO REDUCE THE RISK OF ELECTRIC SHOCK: ALWAYS TURN OFF BEFORE CLEANING OR STORING.

HOW IT WORKS:

Renew mode

uses a medium level of vibrations and red LED light to help reduce the appearance of fine lines and wrinkles.

Restore mode

uses higher vibration levels, blue LED light and heat to help skin appear clearer.

Brighten mode


uses low levels of vibrations and green LED light to help improve the appearance of skin tone.

HOW TO USE:

Device movement

When using the device on your neck, we recommend these movements:

1. From the chin upwards towards your outer cheek.
2. From the upper cheek upwards around your temples.
3. From the bottom of your neck upwards towards your chin.
4. From the bottom of your forehead towards your hairline.

Renew mode

1. Hold down the power button for 1 second to turn on the device and enter the RENEW mode.
2. Apply moisturiser, serum, or lotion to the skin, massage until fully absorbed.
3. Ensure your fingers are touching the silver strip around the handle or the silver ring around the neck to enable the touch sensor.
4. The device will begin vibrating once the tip touches your skin.
5. Glide the massage head along your face and neck for at least 2 minutes.

Restore mode

1. Hold down the power button for 1 second to turn the device on, then quickly press the power button ONCE to enter the RESTORE mode.
2. Apply moisturiser, serum, or lotion to the skin while you wait for the device to heat up. This will take at least 2 minutes. Massage until fully absorbed.
3. Ensure your fingers are touching the silver strip around the handle or the silver ring around the neck to enable the touch sensor.
4. The device will begin vibrating once the tip touches your skin.
5. Glide the massage head along your face and neck for at least 2 minutes.

Brighten mode

1. Hold down the power button for 1 second to turn the device on, then quickly press the power button TWICE to enter the BRIGHTEN mode.
2. Apply moisturiser, serum, or lotion to the skin, massage until fully absorbed.
3. Ensure your fingers are touching the silver strip around the handle or the silver ring around the neck to enable the touch sensor.
4. The device will begin vibrating once the tip touches your skin.
5. Glide the massage head along your face and neck for at least 3 minutes.

SAFETY INFORMATION :

Red Light products are not intended to diagnose, treat, cure, or prevent any disease. Our products are low-risk, general wellness / fitness products that do not require medical clearance. They are designed as "General Wellness" home use products only.

SAFETY CONSIDERATIONS:

Use of a Red Light skincare device in any way other than those described in this manual voids

the warranty and could result in serious injury.

The STYLPRO Fabulous Firmer is an electrical device. To avoid electric shock, adhere to the following safety considerations:

- Do not get your device wet or submerge your device in water.
- Do not wash electrical parts with water or liquids.
- Do not directly touch the device or switches with wet hands.
- Do not apply or touch the device to wet or damp skin.
- Do not pull, carry, or lift the device by its charging cord. If the charging cord is damaged, do not use the device.
- Do not use the device if it has been dropped, left outside, or has encountered water or liquids.
- Do not use on or near heated surfaces.
- Do not modify, or attempt to modify or repair, this equipment. There are no user serviceable parts.
- Do not use or store your device in hot temperatures that are or can exceed 54°C (130°F).

BATTERIES:

THE PRODUCT CONTAINS LITHIUM-ION BATTERY. RISK OF FIRE AND BURNS, DO NOT OPEN, CRUSH, HEAT ABOVE 60°C/140 °F OR INCINERATE. DO NOT SHORT CIRCUIT. IF THE LITHIUM-ION BATTERY BULGES SEVERELY, DISCONTINUE USE. KEEP SMALL CELLS WHICH ARE CONSIDERED SWALLOWABLE OUT OF REACH OF CHILDREN. SWALLOWING MAY LEAD TO BURNS, PERFORATION OF SOFT TISSUE, AND DEATH. SEVERE BURNS CAN OCCUR WITHIN 2 HOURS OF INGESTION. IN CASE OF INGESTION OF A CELL, SEEK MEDICAL ATTENTION PROMPTLY. NEVER RECHARGE NON-RECHARGEABLE BATTERIES. THIS PRODUCT HAS A BUILT-IN RECHARGEABLE BATTERY WHICH CANNOT BE REPLACED BY THE USER. THE REMOVAL OR REPLACEMENT OF THE RECHARGEABLE BATTERY MAY ONLY BE CARRIED OUT BY THE MANUFACTURER OR HIS CUSTOMER SERVICE OR BY A SIMILARLY QUALIFIED PERSON IN ORDER TO AVOID HAZARDS. WHEN DISPOSING OF THE PRODUCT, IT SHOULD BE NOTED THAT THIS PRODUCT CONTAINS A RECHARGEABLE BATTERY.

IMPORTANT SAFETY INSTRUCTIONS:

WHEN USING ELECTRICAL APPLIANCES, BASIC SAFETY PRECAUTIONS SHOULD ALWAYS BE TAKEN, INCLUDING THE FOLLOWING: READ ALL INSTRUCTIONS BEFORE USING APPLIANCE. SAVE THESE INSTRUCTIONS. KEEP AWAY FROM WATER. THIS PRODUCT IS NOT A TOY. DO NOT STARE DIRECTLY INTO THE LIGHTS. USE INDOORS ONLY. FOR DOMESTIC HOUSEHOLD USE ONLY. DO NOT ATTEMPT TO OPEN THE BACK OF THE DEVICE. IF REQUIRED, CONSULT A PROFESSIONAL. KEEP AWAY FROM SOURCES OF DIRECT HEAT AND SUNLIGHT. KEEP AWAY FROM CHILDREN.

WARNINGS: 1. Please turn the device OFF after use. **2.** This product should not be used by; children, pregnant women, patients with epilepsy, tumors and cancer patients, patients with cardiac pacemakers or similar. **3.** Do not place any part of the product near heat sources such as radiators, cookers, etc do not lay on a soft surface such as a bed as device may overheat. **4.** Keep the device away from corrosive liquids. **5.** If you feel abnormalities in your skin during use or after use, please consult your medical specialist. **6.** Do not tamper with any internal components, there are no user serviceable parts inside the unit. Do not attempt to repair the product yourself. **7.** Keep away from water. This product is not a toy. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision. Please store with care and out of reach of children and pets. **8.** To protect against electrical hazards, do not use while bathing or in the shower. **9.** Close supervision is necessary when any appliance is used by or near children or persons made weak or disabled by illness or injury. **10.** Do not use on sensitive or broken skin. Do not use on swollen or sunburnt skin or if you have any other skin condition. **11.** Remove all metals and solid objects, such as earrings, necklaces, bracelets, watches etc before use. Remove, glasses, false eyelashes and contact lenses before use. **12.** Do not operate this appliance if it malfunctions or has been dropped on the floor or in water or damaged in any manner. **13.** Do not directly touch the eyeball. **14.** Do not share this product with anyone else. **15.** Please do not use during charging.

LIMITED TWO-YEAR WARRANTY: OUR GOODS COME WITH GUARANTEES THAT CANNOT BE EXCLUDED UNDER THE EUROPEAN CONSUMER LAW. EU DIRECTIVE 1999/44/EC, THIS DIRECTIVE STATES 'A TWO-YEAR GUARANTEE APPLIES FOR THE SALE OF ALL CONSUMER GOODS EVERYWHERE IN THE EU. YOU ARE ENTITLED TO A REPLACEMENT OR REFUND FOR A MAJOR FAILURE YOU ARE ALSO ENTITLED TO HAVE THE GOODS REPAIRED OR REPLACED IF THE GOODS FAIL TO BE OF ACCEPTABLE QUALITY AND THE FAILURE DOES NOT AMOUNT TO A MAJOR FAILURE. STYLPRO™ WILL REPAIR OR REPLACE (AT OUR OPTION) YOUR APPLIANCE FREE OF CHARGE FOR TWENTY-FOUR MONTHS FROM THE DATE OF PURCHASE, IF THE APPLIANCE IS DEFECTIVE IN WORKMANSHIP OR MATERIALS. THIS WARRANTY IS AN ADDITIONAL BENEFIT AND DOES NOT AFFECT YOUR LEGAL RIGHTS. TO OBTAIN SERVICE UNDER THIS WARRANTY, PLEASE READ THE WARRANTY LEAFLET INCLUDED OR CONTACT STYLPRO™ CUSTOMER SERVICE AT TEAMTOM@STYLIDEAS.COM

+44 (0) 20 3225 5559. TO MAKE A WARRANTY CLAIM YOU MUST OBTAIN YOUR PROOF OF PURCHASE.

IMPORTANT: THIS WARRANTY DOES NOT COVER DAMAGE RESULTING FROM ACCIDENTS, MISUSE, ALTERATIONS, ABUSE, LACK OF REASONABLE CARE OR USE OTHER THAN DESCRIBED IN THESE INSTRUCTIONS, OR PERFORMANCE OF SERVICE OR REPAIRS BY AN UNAUTHORISED SERVICE CENTRE. THIS WARRANTY IS VALID ONLY ON MERCHANDISE SOLD BY AUTHORISED RETAILERS. NOTE, THE ITEM MUST BE THOROUGHLY CLEANED PRIOR TO RETURN.

STORAGE
ENSURE DEVICE IS 'OFF' BEFORE STORAGE. ENSURE THE DEVICE IS UNPLUGGED FROM ANY POWER SUPPLY. PLACE IN COOL, DRY AREA OUT OF REACH OF CHILDREN, WHERE IT WILL NOT BE DAMAGED.

MAINTENANCE
YOUR STYLPRO™ Fabulous Firmer IS VIRTUALLY MAINTENANCE FREE. TO CLEAN TURN THE APPLIANCE OFF AND Wipe EXTERIOR WITH A DAMP CLOTH. IF ANY ABNORMAL CONDITION OCCURS, TURN OFF & CONTACT STYLPRO™ CUSTOMER SERVICES ON +44 (0) 20 3225 5559.

CORRECT DISPOSAL OF THIS PRODUCT: WHEN THE PRODUCT HAS REACHED THE END OF ITS USAGE, PLEASE DO NOT DISPOSE WITH REGULAR HOUSEHOLD WASTE. RECYCLING IS THE SAFEST AND MOST RESPONSIBLE WAY TO DISPOSE OF ELECTRICAL ITEMS THAT COME TO THE END OF THEIR LIFE CYCLE. PLEASE DISPOSE OF IT ACCORDING TO YOUR LOCAL REGULATIONS. THE WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT (WEEE) DIRECTIVE IN EUROPEAN COUNTRIES STIPULATES CORRECT RECYCLING PRACTICES. DO NOT DISPOSE IN REGULAR HOUSEHOLD WASTE.

CE MARKING: UK/EU - DECLARATION OF CONFORMITY: STYLIDEAS LTD DECLARES THAT STYLPRO COMPLIES WITH THE ESSENTIAL REQUIREMENTS AND OTHER RELEVANT PROVISIONS. A COPY OF THE DECLARATION OF CONFORMITY IS AVAILABLE ON REQUEST. STYLIDEAS LTD, AMSHOLD HOUSE, GOLDINGS HILL, LOUGHTON, ESSEX, IG10 2RW, UK. OBELIS S.A, BD. GENERAL WAHIS 53, B-7030 BRUSSELS, BELGIUM.

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